

HACH-LANGE Electrochemistry

Determination of Chloride in alcohols using an ISE (Ion selective electrode)

Application APP-ECH-0022
DOC042.52.20138.May13

Introduction


Bio-alcohols are used in bio-fuels. The composition of those fuels has to be controlled to avoid ingredients, which may corrode the metal housing of an engine. Chloride (Cl^-) is a very corrosive ion and must be absent, because it can initiate a corrosion process, what finally leads to damage.

The Hach HQ40d Multimeter with ISECL18101 ion selective probe for Cl^- ions, is a perfect system to analyze for Cl^- in alcohols. Normally, electrochemical sensors are designed for aqueous solutions only and organic solutions like alcohol need a specific reference system. The IntelliCal ISECL18101 works without any modification directly in the alcoholic sample.

Calibration and sample measurement

Prepare at least 4 standard solutions of an alcohol / water mixture what is identical to the composition of the sample. For Methanol a 100% Methanol solution can be used to dissolve Sodium Chloride (NaCl). Ethanol needs some water to dissolve NaCl , because pure Ethanol will not dissolve NaCl . Therefore a 95% Ethanol + 5% DI water mixture can be used for the calibration.

Start with a stock solution of 1000 mg/L Cl^- , what is 1.655 g NaCl salt dissolved in 1 L Methanol. Take corresponding volumes to dilute the stock solution to the required lower Cl^- standards. The table on the right shows a measurement of the ISECL18101 in 6 standards. Using 1, 10 and 100 mg/L for calibration the calculated slope is -48.5 mV/pX and offset +107mV. The ISECL18101 reacts linear over the complete range between 1 and 100 mg/l Cl^- in Methanol.




| Cl ⁻ conc. (mg/l) | Log (c) | mV measured |
|---------------------------------|----------|-------------|
| 1 | 0,00 | 107,8 |
| 2 | 0,30 | 94,1 |
| 5 | 0,70 | 72,0 |
| 10 | 1,00 | 57,7 |
| 20 | 1,30 | 44,0 |
| 100 | 2,00 | 11,6 |

The Chloride analysis in Ethanol is different to Methanol, because Ethanol is less polar than Methanol or Water. Below the table shows the measurements of ISECL181 in a 95 Ethanol solution with different amount of Cl^- .

A stock solution of 1,655 g NaCl / 1 L was prepared with 95% Ethanol + 5 % DI water. Compared to Methanol the calibration curve of Cl^- in Ethanol is not linear, but a cubic function. Therefore the sensitivity at lower Cl^- concentrations decreases faster than in Methanol.

The slope calculated from the 1, 10, 100 mg/L Cl^- standards is -43.9 mV/pX and offset 72.5 mV.

The ISECL18101 reacts almost linear over the range between 1 and 100 mg/l Cl^- in Ethanol.



| Cl- concentration (mg/l) | LOG (C) | mV |
|--------------------------|-----------|-------|
| 0,1 | -1 | 92,4 |
| 1 | 0 | 71,0 |
| 10 | 1 | 31,5 |
| 100 | 2 | -16,9 |

Maintenance

After each measurement rinse the ISE thoroughly with DI water and store the probe in the storage solution.