



Solución estándar de cloruro de sodio, 10 000 $\mu\text{S}/\text{cm}$, 500 mL

# de producto:	2972249
USD Precio:	Contacto Hach
Disponible	

Conductivity standard solutions for laboratory and field analyses

Para comprobar la exactitud de las mediciones de conductividad. 10 000 $\mu\text{S}/\text{cm}$. Botella de 500 mL.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate conductivity measurements is good care and calibration of the conductivity cell.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Especificaciones

Cantidad:	500 mL
Concentración:	5,556 g/L $\pm 0,03$ g/L NaCl
Condiciones de almacenamiento:	10 - 25 °C (protect from light)
Descripción:	Cloruro sódico, 10 000 $\mu\text{S}/\text{cm}$, 500 ml
Parámetro:	Conductividad
Plataforma:	Solución
Tipo de envase:	Bottle