



Solución tampón, pH 7,00, cert., 1 L

# de producto:	LZW9467.97
USD Precio:	Contacto Hach
Disponible	

pH Buffer standard solutions for laboratory and field analyses

Solución estándar tampón de pH con certificado de análisis. Solución de pH 7,00. Botella de 1 L.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Especificaciones

Cantidad:	1 L
Condiciones de almacenamiento:	10 - 25 °C (protect from light)
Descripción:	pH Buffer Calibration reagent
Notas especiales:	Observed shelf-life after opening, following recommended storage and use: 6 months
Tipo de envase:	Bottle
Trazable a NIST:	Sí
Valor(es) pH:	7,00
¿Codificada por colores?:	No